

Methamphetamines

Waters Corporation



For forensic toxicology use only.

This is an Application Brief and does not contain a detailed Experimental section.


Abstract

This application brief highlights the analysis of methamphetamines using XTerra Columns.


Introduction

The compounds analyzed in this study are:

1. 3,4-methylenedioxyamphetamine (MDA)
2. 3,4-methylenedioxymethamphetamine (MDMA)



3,4-Methylenedioxyamphetamine (MDA)




3,4-Methylenedioxymethamphetamine (MDMA)

Experimental

Conditions

Column:	XTerra RP ₁₈ 3.0 x 100 mm, 3.5 μm
Part number:	186000420
Mobile phase A:	5 mM NH ₄ HCO ₃ , pH 9.6
Mobile phase B:	ACN
Flow rate:	0.6 mL/min
Isocratic mobile phase composition:	80%A; 20%B
Injection volume:	5 μL
Sample concentration:	40 pg/μL
Temperature:	40 °C
Detection:	Fluorescence 280ex 326em Gain 100 TC = 1s
Instrument:	Alliance 2690, 2475
Compounds	USP Tailing
3,4-methylenedioxyamphetamine (MDA)	0.95
3,4-methylenedioxymethamphetamine (MDMA)	1.10

Results and Discussion



Featured Products

- [Alliance HPLC System <https://www.waters.com/534293>](https://www.waters.com/534293)
- [2475 Fluorescence \(FLR\) Detector <https://www.waters.com/514434>](https://www.waters.com/514434)

WA20738.065, June 2002



© 2021 Waters Corporation. All Rights Reserved.