

## LC-MS Analysis of Vitamin B<sub>12</sub>

---

Hongjian Jiang, Zhao Gui Pin, Kate Yu

Waters Corporation

---

### Abstract

This application note describes a simple LC-MS method for the analysis of a Vitamin B<sub>12</sub> standard.

### Benefits

With the use of the Waters Atlantis C<sub>18</sub> Column, Vitamin B<sub>12</sub> was well retained, with a good peak shape

---

### Introduction

Vitamin B<sub>12</sub> (cyanocobalamine) is a water-soluble vitamin that has important physiological function in the human body. Vitamin B<sub>12</sub> has the ability to make folic acid available to bone marrow, which appears to be necessary for red blood cell formation. A Vitamin B<sub>12</sub> deficiency leads to a degeneration of both the sensory and motor columns in the spinal cord with loss of sensation and paralysis. As a result, people need approximately 1–2 µg of Vitamin B<sub>12</sub> per day.

The determination of Vitamin B<sub>12</sub> is important, but can be challenging largely because of its chemical instability and the complexity of the matrices in which it is usually found. Ideally, methods for Vitamin B<sub>12</sub> analysis should be simple, selective and sensitive to overcome the above issues. The most commonly reported method for Vitamin B<sub>12</sub> analysis is reversed-phase HPLC separation with various detection schemes such as

---

UV, fluorescence, chemiluminescence, etc. However, most of these methods are tedious and time consuming. Phosphate buffers are usually called for in order to obtain sufficient LC separation.

In this application note, we describe a simple LC-MS method for the analysis of a Vitamin B<sub>12</sub> standard. A Waters Atlantis C<sub>18</sub> Column, with superior retention for polar compounds and the ability to operate in 100% aqueous conditions, was used for the separation. A Waters Micromass ZQ 4000 single quadrupole mass spectrometer was used for detection. With the selectivity and sensitivity offered by the MS detector, the method simply used a binary acetonitrile/water gradient without the need for a buffer or ion pairing reagents.

---

## Experimental

### LC Conditions

LC system:	Alliance HT Separations Module
Column:	Atlantis C <sub>18</sub> 2.1 x 150 mm, 3.5 µm
Flow rate:	0.2 mL/min
Mobile phase:	Acetonitrile (A) Water (B)
Injection volume:	10 µL

### Gradient

Time(min)	%A	%B	Curve
0	0	100	1
1	0	100	1
8	40	60	6

Time(min)	%A	%B	Curve
10	60	40	6
11	0	100	11

## MS Conditions

Mass spectrometer:	Micromass ZQ 4000 Mass Spectrometer
Ion mode:	ESI+
Capillary voltage:	3.5 kV
Cone voltage:	80 V
Source temp.:	105 °C
Desolvation temp.:	180 °C
Desolvation gas flow:	365 L/Hour
Cone gas flow:	50 L/Hour
Inter channel delay:	0.02 s
Inter scan delay:	0.02 s
Dwell time:	0.08 s

---

## Results and Discussion

Full scan spectra of Vitamin B<sub>12</sub> were obtained during the optimization of the MS parameters, as shown in Figure

---

1. Acquisition parameters were optimized in order to select the proper ion for selected ion monitoring (SIR) experiments and to determine an optimum cone voltage. The protonated molecule  $[M+H]^+$  at  $m/z$  1356 was obtained at a cone voltage of 50 V with good intensity ( $4.41e5$ ). However, the Vitamin B<sub>12</sub> showed fragmentation into many other ions, which reduced the intensity of the  $[M+H]^+$  ion. The spectrum at 80 V was cleaner with higher ion intensity ( $6.14e5$ ) at  $m/z$  913. As a result, the ion at  $m/z$  913 was chosen to be the ion to monitor for the SIR experiments and cone voltage of 80 V was chosen as the optimum for the ion at  $m/z$  913.




Figure 1. Full scan MS spectra for Vitamin B<sub>12</sub> at different cone voltages.

Figure 2 shows the calibration curve of Vitamin B<sub>12</sub>. The linear range was from 5 to 4000 ppb. With the Waters Atlantis C<sub>18</sub> Column, this water-soluble vitamin was well retained with an excellent peak shape. Figure 3 shows the SIR chromatograms of Vitamin B<sub>12</sub> at *m/z* 913 at different concentrations. The signal-to-noise ratio at 2 ppb was 5, and at 5 ppb was 50.




Figure 2. Calibration curve of Vitamin B<sub>12</sub>.




Figure 3. SIR chromatograms of Vitamin B<sub>12</sub> at m/z 913.

---

## Conclusion

A simple LC-MS quantification method for Vitamin B<sub>12</sub> was described. With the use of the Waters Atlantis C<sub>18</sub> Column, Vitamin B<sub>12</sub> was well retained, with a good peak shape. A simple binary gradient was used without adding buffer. The quantification linear range was from 5n to 4000 ppb. The limit of detection with a 10 μL injection was 2 ppb.

---

## Featured Products

Alliance HPLC <<https://www.waters.com/514248>>

720000758, September 2003

© 2022 Waters Corporation. All Rights Reserved.