

Application Note

## Aflatoxins in Peanuts

---

YU Yanbin, WAN Shuwei, TAN Peigong, WANG Xiaogang, MIAO Zaijing

Ministry of Agriculture, Qingdao, Environmental Protection Agency, Qingdao, Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

---

### Abstract

This application brief describes the analysis of aflatoxins in peanuts.

---

## Introduction

This application brief describes the analysis of aflatoxins in peanuts.

---


## Experimental

### Pretreatment

1. Add 5 g of sodium hydroxide to 20 g of homogenized sample, followed by 30 mL of n-hexane.
2. Add 100 mL 60% aqueous methanol and homogenize.
3. Ultrasonicate for 30 minutes.
4. Filter sample through 15 cm filter paper.
5. Take 1 mL aliquot from 60% methanol layer for SPE cleanup.

## SPE Procedure

Oasis® HLB 1cc/30 mg




## LC Conditions

|                       |  |
|-----------------------|--|
| Instrument:           | Alliance HPLC 2695 System                |
| Column:               | Symmetry Shield RP18, 4.6 x 150 mm, 5 µm |
| Flow rate of iodine:  | 0.2 mL                                   |
| Flow rate:            | 1 mL/min                                 |
| Mobile phase:         | A. methanol<br>B. water                  |
| Isocratic gradient:   | 35% A: 65% B, for 20 minutes             |
| Column temperature:   | 30 °C                                    |
| Derivatization temp.: | 80 °C                                    |

|                                     |   |
|-------------------------------------|---|
| Excitation wavelength:              | 365 nm  |
| Emission wavelength:                | 455 nm  |
| Post-column derivatization reagent: | Dissolve 200 mg iodine in 10 mL methanol, top up 1000 mL with water |
| Detector:                           | 2475 Multi Wavelength Fluorescence                                  |

---

## Results and Discussion



Matrix interference is greatly reduced when sample is cleaned up by using Oasis HLB SPE cartridge.

| Analyte      | Recovery % | Detection (p/μg kg) |
|--------------|------------|---------------------|
| Aflatoxin G2 | 101± 7.18  | 0.11                |
| Aflatoxin G1 | 72.8±3.63  | 0.20                |
| Aflatoxin B2 | 97.5±5.48  | 0.12                |
| Aflatoxin B1 | 68.8±5.48  | 0.24                |

Results of B1, B2, G1, G2 in peanuts (n=5)

---

## Featured Products

Alliance HPLC System <<https://www.waters.com/534293>>

720002584, April 2008

© 2021 Waters Corporation. All Rights Reserved.