


## Gradient Separation of Histidine Dipeptides on ACQUITY UPLC BEH HILIC

---

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

---

### Abstract


This application brief demonstrates the gradient separation of histidine dipeptides on ACQUITY UPLC BEH HILIC Columns.

---


## Introduction

The compounds used in this study are:


1. Creatinine
2. Creatine
3. Anserine
4. Carnosine




Creatinine



Creatine



Anserine



Carnosine

---

## Experimental

### Test Conditions

Column:

ACQUITY UPLC BEH HILIC, 2.1 x 50 mm, 1.7 μm

Part Number:

186003460

Mobile Phase A:	50/50 ACN/10 mM ammonium formate, w/ 0.125% HCOOH, pH 3.0
Mobile Phase B:	95/5 ACN/10 mM ammonium formate, w/ 0.125% HCOOH, pH 3.0
Flow Rate:	0.5 mL/min
Injection Volume:	5.0 µL
Sample Diluent:	75:25 ACN:MeOH
Sample Concentration:	Creatinine 1µg/mL; Carnosine 5 µg/mL; Anserine 5 µg/mL; Creatine 5 µg/mL
Column Temperature:	30 °C
Weak Needle Wash:	ACN/H <sub>2</sub> O 95/5
Instrument:	Waters ACQUITY UPLC with ACQUITY SQD

## Gradient

Time (min)	Profile
	%A
0.00	0.1
5.00	99.9
5.01	0.1
6.00	0.1

## MS Conditions

Ionization Mode: ES+

Capillary: 2.5 kV

Cone: 20 V (Carnosine; Creatinine, Anserine);  
25 V (Creatine)

Source Temperature: 120 °C

Desolvation Temperature: 400 °C

Desolvation Gas Flow: 800 L/Hr


Cone gas Flow: 5 L/Hr

SIR *m/z*: 227.1 *m/z* (Carnosine); 132.1 *m/z* (Creatine);  
114.05 *m/z* (Creatinine); 241.1 *m/z* (Anserine)

Dwell Time: 0.1 s

---

## Results and Discussion



---

## Featured Products

ACQUITY UPLC System <<https://www.waters.com/514207>>

WA60139, August 2009

© 2021 Waters Corporation. All Rights Reserved.