


## Gradient Separation of Histidine Dipeptides on ACQUITY UPLC BEH HILIC

---

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

### Abstract

This application brief demonstrates the gradient separation of histidine dipeptides on ACQUITY UPLC BEH


---

HILIC Columns.


## Introduction

The compounds used in this study are:


1. Creatinine
2. Creatine
3. Anserine
4. Carnosine




Creatinine



Creatine



Anserine



Carnosine

## Experimental

### Test Conditions

---

Column: ACQUITY UPLC BEH HILIC, 2.1 x 50 mm, 1.7  $\mu$ m

Part Number: 186003460

Mobile Phase A: 50/50 ACN/10 mM ammonium formate, w/ 0.125% HCOOH, pH 3.0

Mobile Phase B: 95/5 ACN/10 mM ammonium formate, w/ 0.125% HCOOH, pH 3.0

Flow Rate: 0.5 mL/min

Injection Volume: 5.0  $\mu$ L

Sample Diluent: 75:25 ACN:MeOH

Sample Concentration: Creatinine 1 $\mu$ g/mL; Carnosine 5  $\mu$ g/mL; Anserine 5 g/mL; Creatine 5  $\mu$ g/mL

Column Temperature: 30  $^{\circ}$ C

Weak Needle Wash: ACN/H<sub>2</sub>O 95/5

Instrument: Waters ACQUITY UPLC with ACQUITY SQD

## Gradient

Time (min)	Profile
	%A
0.00	0.1
5.00	99.9


---

Time (min)	Profile
5.01	0.1
6.00	0.1

## MS Conditions

Ionization Mode:	ES+
Capillary:	2.5 kV
Cone:	20 V (Carnosine; Creatinine, Anserine); 25 V (Creatine)
Source Temperature:	120 °C
Desolvation Temperature:	400 °C
Desolvation Gas Flow:	800 L/Hr
Cone gas Flow:	5 L/Hr
SIR <i>m/z</i> :	227.1 <i>m/z</i> (Carnosine); 132.1 <i>m/z</i> (Creatine); 114.05 <i>m/z</i> (Creatinine); 241.1 <i>m/z</i> (Anserine)
Dwell Time:	0.1 s

## Results and Discussion



## Featured Products

· [ACQUITY UPLC System <https://www.waters.com/514207>](https://www.waters.com/514207)

WA60139, August 2009



© 2021 Waters Corporation. All Rights Reserved.