

Analysis of Morphine and Its Polar Metabolites Using Atlantis T3 Columns

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

Abstract

This application brief demonstrates analysis of morphine and its polar metabolites using Atlantis T3 Columns.

Introduction

The compounds used in this study are -

1. 10-Hydroxymorphine
2. Morphine-3 β -D-glucuronide
3. Morphine-6 β -D-glucuronide
4. Morphine
5. Morphine N-oxide
6. 6-Acetylmorphine

Experimental

Test Conditions

Column:	Atlantis T3, 2.1 x 50 mm, 3 μ m
Part number:	186003717
Mobile phase A:	0.1% FA in water
Mobile phase B:	ACN
Flow rate:	0.5 mL/min
Injection volume:	15 μ L
Column temperature:	45 $^{\circ}$ C


MRM: morphine 286 > 201
morphine-3 β -D-glucuronide 462 > 286
morphine-6 β -D-glucuronide 462 > 286
morphine N-oxide 302 > 162
6-acetylmorphine 328 > 165
10-hydroxymorphine 302 > 58

Instrument: Waters ACQUITY UPLC System

Gradient

Time (min)	Profile	
	%A	%B
0.00	98	2
5.00	2	98

Results and Discussion



Featured Products

· [ACQUITY UPLC System <https://www.waters.com/514207>](https://www.waters.com/514207)

WA60210, January 2009



© 2021 Waters Corporation. All Rights Reserved.