

## Gradient Separation of Morphine and Metabolites on XBridge HILIC

---

Waters Corporation



This is an Application Brief and does not contain a detailed Experimental section.

---

### Abstract


This application brief demonstrates the gradient separation of morphine and metabolites on Xbridge HILIC Columns.

---


## Introduction

The compounds used in this study are:


1. 6-Acetylmorphine
2. Morphine
3. Morphine-3 $\beta$ -D-glucuronide



6-Acetylmorphine



Morphine



Morphine-3 $\beta$ -D-glucuronide

---

## Experimental

## Method Conditions

Column:	XBridge HILIC, 2.1 x 50 mm, 3.5 µm
Part number:	186004432
Mobile phase A:	10 mM NH <sub>4</sub> COOH in H <sub>2</sub> O, 0.125% HCOOH in 50:50 ACN:H <sub>2</sub> O
Mobile phase B:	10 mM NH <sub>4</sub> COOH in H <sub>2</sub> O, 0.125% HCOOH in 90:10 ACN:H <sub>2</sub> O
Flow rate:	0.6 mL/min
Injection volume:	5 µL
Sample concentration:	25 ng/mL each
Sample diluent:	75:25 ACN:MeOH with 0.2% HCOOH
Column temperature:	30 °C
Detection:	UV @ 280 nm
Sampling rate:	20 points/sec
Time constant:	0.1
Instrument:	Waters ACQUITY UPLC with ACQUITY PDA


## Gradient:

Time(min)	Profile
	%A
0.00	0.1

Time(min)	Profile
1.05	0.1
4.35	99.9
4.50	0.1
6.00	0.1

---

## Results and Discussion



---

## Featured Products

ACQUITY UPLC System <<https://www.waters.com/514207>>

ACQUITY UPLC PDA Detector <<https://www.waters.com/514225>>

WA64079, August 2009

© 2021 Waters Corporation. All Rights Reserved.