

Note d'application

## Quantitative Analysis of Barbiturates in Urine Using UPLC-MS/MS

---

Robert Lee, Allan Traynor, Jane LeCount, Michelle Wood

Waters Corporation, Concateno plc



*For forensic toxicology use only.*

---

## Abstract

Analyzing barbiturates by UPLC-MS/MS, the use of a very simple sample dilution step eliminates both the liquid-liquid extraction and post-extraction derivatization steps that are required for GC-MS analysis. The elimination of the extraction step would reduce the time taken to prepare a typical batch of samples by more than 50%.

### Benefits

- Elimination of extraction step
- Elimination of derivatization step prior to analysis
- Improved sample throughput

---

## Introduction

Barbiturates act as central nervous system depressants producing effects ranging from mild sedation to general anesthesia. They have largely been replaced by benzodiazepines as prescription medicines, owing to their relatively low therapeutic index and their high potential for dependence. However, it is known that the use of barbiturates is still common in certain regions of Eastern Europe;<sup>1</sup> consequently, their analysis is still of key importance in both forensic analysis and workplace drug testing.

Barbiturates have traditionally been measured by GC.<sup>2,3</sup> The arrival of newer technologies into the modern laboratory, such as UPLC-MS/MS, often leads to an overall requirement to consolidate analytical methods and transfer existing methodologies to the newer platforms. Furthermore, UPLC-MS/MS permits the development of more sensitive techniques.

We report a quantitative method based on simple dilution and UPLC-MS/MS. The method has been verified, and its performance evaluated using authentic samples. Data were compared to results obtained with a traditional method that used liquid-liquid extraction followed by derivatization and analysis by GC-MS.

---

## Experimental

### Sample description

Phenobarbital and thiopental were purchased from Sigma Aldrich (Dorset, UK) and dissolved in methanol to 1 mg/mL. All other barbiturates (1 mg/mL) and deuterated internal standards (ISTDs) at 0.1 mg/mL were obtained as certified standard solutions from LGC Standards (Teddington, UK). Deuterated internal standards were not available for all of the barbiturates.

Quality control reference urine samples (Bio-Rad Liquichek Urine Toxicology Control: C2, C3, C4, S1, and S2) were obtained from Bio-Rad Laboratories Ltd (Hemel Hempstead, UK).

Urine samples for method development were obtained from donors at Waters Corporation.

Nineteen samples containing pre-analyzed barbiturates were obtained from Concateno, London, UK.

### Sample preparation

Urine, either sample or calibrator, was centrifuged at 13,000 rpm for 5 min, then 50 µL was transferred to a Waters maximum recovery vial and diluted with 950 µL water containing 25 ng of each available ISTD.

### UPLC Conditions

|                   |  |
|-------------------|--|
| System:           | ACQUITY UPLC   |
| Column:           | ACQUITY BEH C <sub>18</sub> 2.1 x 100 mm with BEH C <sub>18</sub> 2.1 x 5 mm Vanguard pre-column |
| Column temp.:     | 50 °C  |
| Sample temp.:     | 5 °C   |
| Injection volume: | 15 µL (PLNO)   |
| Strong wash:      | 0.001% formic acid in acetonitrile   |
| Weak wash:        | 0.001% formic acid in water  |

|                 |  |
|-----------------|--|
| Flow rate:      | 400 $\mu$ L/min  |
| Mobile phase A: | 0.001% formic acid in water  |
| Mobile phase B: | 0.001% formic acid in acetonitrile   |
| Gradient:       | Hold at 5% B for 0.5 min, then switch to 27.5% B, hold until 4 min, then switch to 35% B, hold until 5.25 min, then switch to 90% B, hold until 6.25 min, then switch to 5% B. |

## MS conditions

|                    |  |
|--------------------|--|
| Mass Spectrometer: | Xevo TQ  |
| Ionization mode:   | ESI negative   |
| Capillary voltage: | 2.75 kV  |
| Cone voltage:      | 25 V   |
| Collision energy:  | 12 eV  |
| Desolvation temp.: | 500 °C   |
| Desolvation gas:   | 1000 L/h   |
| Cone gas:          | 25 L/h   |
| Acquisition mode:  | Multiple reaction monitoring (MRM), as shown in Table 1. |

## Data management

MassLynx v4.1 incorporating TargetLynx Application Manager

| Barbiturate      | Precursor ion<br>( <i>m/z</i> ) | Product ion<br>( <i>m/z</i> ) | ISTD             |
|------------------|---------------------------------|-------------------------------|------------------|
| Amobarbital      | 225.1                           | 182.0                         | Pentobarbital-d5 |
| Butabarbital     | 211.1                           | 168.0                         | Phenobarbital-d5 |
| Butalbital       | 223.0                           | 180.0                         | Butalbital-d5    |
| Pentobarbital    | 225.1                           | 182.0                         | Pentobarbital-d5 |
| Phenobarbital    | 231.1                           | 188.0                         | Phenobarbital-d5 |
| Secobarbital     | 237.1                           | 194.1                         | Secobarbital-d5  |
| Thiopental       | 241.1                           | 57.9                          | Secobarbital-d5  |
| Butalbital-d5    | 228.0                           | 185.0                         |                  |
| Pentobarbital-d5 | 230.1                           | 186.9                         |                  |
| Phenobarbital-d5 | 236.1                           | 193.0                         |                  |
| Secobarbital-d5  | 242.1                           | 199.1                         |                  |

*Table 1. MRM transitions for analytes and ISTDs.*

## Results and Discussion

### Method verification

The MRM transitions for all of the barbiturates and ISTDs are shown in Table 1. All were monitored using a single transition. Figure 1 shows a chromatogram of a 500 ng/mL barbiturate-spiked urine.




Figure 1. Chromatogram showing barbiturates spiked into urine at 500 ng/mL.

To investigate linearity for all barbiturates, spiked urine calibrators were prepared at 0, 25, 50, 100, 250, 500, Quantification was performed by integrating the area under the peak for each analyte MRM trace, and refer

2. Pocci R *et al. J Anal Toxicol.* 1992; 16(1): 45-72.

3. Johnsonn LL and Garg U. *Methods Mol Biol.* 2010; 603: 65-74.

## Acknowledgments

Concateno plc, London, UK for supplying the anonymized authentic urine samples.

*A full validation by the user would be necessary prior to adoption in a laboratory.*

---

## Featured Products

ACQUITY UPLC System <<https://www.waters.com/514207>>

TargetLynx <<https://www.waters.com/513791>>

MassLynx MS Software <<https://www.waters.com/513662>>

720004466, September 2014