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
1. Periat A, Boccard J, Veuthey JL, Rudaz S, Guillarme D. Systematic Comparison of Sensitivity between Hydrophilic Interaction Liquid Chromatography and Reversed Phase Liquid Chromatography. *J. Chrom. A*, 2013, 1312 49–57.
2. HILIC Separations, A Practical Guide to HILIC Mechanisms, Method Development and Troubleshooting. Thermo Scientific Technical Guide. <https://assets.thermofisher.com/TFS-Assets/CMD/brochures/TG-21003-HILIC-Separations-TG21003-EN.pdf> Accessed 9-July-2020 <<https://assets.thermofisher.com/TFS-Assets/CMD/brochures/TG-21003-HILIC-Separations-TG21003-EN.pdf%20Accessed%209-July-2020>>
3. Hung, I. F.-N. Triple Combination of Interferon beta-1b, lopinavir-Ritonavir, and Ribavirin in the Treatment of Patients Admitted to Hospital with COVID-19: An Open-label, Randomized, Phase 2 Trial, *Lancet*, 2020, 395, 1695–1704.
4. cLogD值参考网站: <https://chemicalize.com/welcome> <<https://chemicalize.com/welcome>>

5. McCalley D. Understanding and Manipulating the Separation in Hydrophilic Interaction Liquid Chromatography.*J.Chrom.A*, 2017, 1523, 49–71.
 6. Guo Y. Recent Progress in the Fundamental Understanding of Hydrophilic Interaction Chromatography (HILIC).*Analyst*, 2015, 140, 6452-6466.
 7. Periat A, Debrus B, Rudaz S, Guillarme D. Screening of the Most Relevant Parameters for Method Development in Ultra-High Performance Hydrophilic Interaction Chromatography.*J. Chrom.A*, 2013, 1282, 72-83.
 8. Danso D, Langman LJ, Snozek C. LC-MS/MS Quantitation of Ribavirin in Serum and Identification of Endogenous Isobaric Interferences.*Clin.Chim.Acta*.2011, 412, 2332–2335.
-

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