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1. SV Bhat, BA Nagasampagi, M Sivakumar.Chemistry of Natural Products.2005, ISBN 3-540-40669-7.
2. H Koltai, P Poulin, D Namdar.Promoting Cannabis Products to Pharmaceutical Drugs. European J. Pharm.Sci. 2019, 132, 118–120.
3. AA Izzo, F Borrelli, R Capasso, V Di Marzo, R Mechoulam.Non-Psychotropic Plant Cannabinoids: New Therapeutic Opportunities from an Ancient Herb.Trends in Pharmaco.Sci. 2009, 30, (10), 515–527.
4. CM Andre, JF Hausman, G Guerriero.Cannabis sativa: The Plant of the Thousand and One Molecules.Front.in Plant Sci.2016, 7(19).
5. MA ElSohly, D Slade.Chemical Constituents of Marijuana: The Complex Mixture of Natural Cannabinoids.Life Sciences. 2005, 78, 539–548.
6. R Brenneisen.Chemistry and Analysis of Phytocannabinoids and Other Cannabis Constituents. Marijuana and the Cannabinoids, 2007, 17–49.doi:10.1007/978-1-59259-947-9_2.
7. C Citti, P Linciano, F Russo, L Luongo, M Iannotta, S Maione, A Lagana, AL Capriotti, F Forni, MA Vandelli, G Gigli, G Cannazza.A Novel Phytocannabinoid Isolated from Cannabis sativa L.with an in vivo Cannabimimetic Activity Higher than Δ^9 -tetrahydrocannabinol: Δ^9 -tetrahydrocannabiphorol. Scientific Reports.2019, 9:20335.

8. MM Radwan, MA ElSohly, D Slade, SA Ahmed, IA Khan, SA Ross, Biologically Active Cannabinoids from High Potency Cannabis sativa. *J. Nat. Prod.* 2009, 906-911.
9. C Citti, D Braghiroli, MA Vandelli, G Cannazza. Pharmaceutical and Biomedical Analysis of Cannabinoids: A Critical Review. *J. Pharm. Biomed. Sci.* 2018, 147, 565-579.
10. MM Lewis, Y Yang, E Wasilewski, HA Clarke, LP Kotra. Chemical Profiling of Medical Cannabis Extracts. *ACS Omega*. 2017, 2, 6091-6103.
11. JT Fishedick, A Hazekamp, T Erkelens, YH Choi, R Verpoorte. Metabolic Fingerprinting of Cannabis sativa L., Cannabinoids and Terpenoids for Chemotaxonomic and Drug Standardization Purpose. *Phytochemistry*. 2010, 71, 2058-2073.
12. A Hazekamp, JT Fishedick. Cannabis – From Cultivar to Chemovar. *Drug Test. Analysis*. 2012, DOI 10.1002/dta.407.
13. A Hazekamp, K Tejkalova, S Papadimitriou. Cannabis: From Cultivar to Chemovar II-A metabolomics Approach to Cannabis Classification. *Cannabis and Cannabinoid Research*. 2016, 1.1, 202–215.
14. P Berman, K Futoran, GM Lewitus, D Mukha, M Benami, T Shlomi, D Meiri. A New ESI-LC/MS Approach for Comprehensive Metabolic Profiling of Phytocannabinoids in Cannabis. *Scientific Reports*. 2018, 8:14280.
15. EB Russo. Taming THC: Potential Cannabis Synergy and Phytocannabinoid-Terpenoid Entourage Effects. *British Journal of Pharmacology*. 2011, 163, 1344-1364.
16. MM Delgado-Povedano, C Sanchez-Carnerero, F Priego-Capote. Untargeted Characterisation of Extracts from Cannabis sativa L., Cultivars by Gas and Liquid Chromatography Coupled to Mass Spectrometry in High Resolution Mode. *Talanta*. 2020, 208, 120384.
17. S Elzinga, J Fishedick, R Podkolinski, JC Raber. Cannabinoids and Terpenes as Chemotaxonomic Markers in Cannabis. *Nat. Prod. Chem. Res.* 2015, 3(4).
18. KW Hillig, A Chemotaxonomic Analysis of Terpenoid Variation in Cannabis. *Biochemical Systematics and Ecology*. 2004, 875-891.
19. JT Fishedick. Identification of Terpenoid Chemotypes Among High (-)-trans- Δ^9 -tetrahydrocannabinol-Producing Cannabis sativa L. Cultivars. *Cannabis and Cannabinoid Research*. 2017, 2.1, 34–47.

20. C Orser, S Johnson, MD Speck, A Hilyard, Ini Afia. Terpenoid Chemoprofiles Distinguish Drug-Type Cannabis sativa L. cultivars in Nevada. Nat. Prod. Chem. Res. 2018, 6(304).
21. J Booth, J Bohlman. Terpenes in Cannabis sativa – From Plant Genome to Humans. Plant Science: An International Journal of Experimental Plant Biology. 2019 Vol: 284, Page: 67–72.
22. G Magagnini, G Grassi, S Kotiranta. The Effect of Light Spectrum on the Morphology and Cannabinoid Content of Cannabis sativa L., Med Cannabis Cannabinoids. 2018, 1:19–27.
23. EA Ibrahim, Wang M, MM Radwan, AS Wanas, CG Majumdar, B Avula, YH Wang, IA. Khan, S Chandra, H Lata, GM. Hadad, RA Abdel Salam, AK. Ibrahim, SA Ahmed, MA. ElSohly. Analysis of Terpenes in Cannabis sativa L. Using GC-MS: Method Development, Validation, and Application. Planta Med. 2019, 85, 431-438.
24. L Eriksson, E Johansson, N Kettaneh-Wold, J Tyrgg, C Wikstrom, S Wold. 2006, Multi- and Megavariate Data Analysis: Part 1 Basic Principles and Applications (2nd ed).
25. G Navarro, K Varani, I Reyes-Resina, et al. Cannabigerol Action at Cannabinoid CB1 and CB2 Receptors and at CB1-CB2 Heteroreceptor Complexes. Front Pharmacol. 2018, 9:632. doi:10.3389/s41370-018-0094-1.
26. Da Cheng Hao, Xiao-Jie Gu, Pei Gen Xiao. 2015, Medicinal Plants: Chemistry, Biology and Omics. Phytochemical and Biological Research of Cannabis Pharmaceutical Resources Chapter 11. 431–464. <http://dx.doi.org/10.1016/B978-0-08-100085-4.00011-6> <<http://dx.doi.org/10.1016/B978-0-08-100085-4.00011-6>> .
27. JD Storey, R Tibshirani, Statistical Significance of Genome Wide Studies. Statistical Significance for Genomewide Studies. Proc Natl Acad Sci U S A. 2003; 100(16):9440–9445. doi:10.1073/ pnas.1530509100.
28. MA ElSohly, MM Radwan, W Gul, S Chandra, A Galal. Phytochemistry of Cannabis sativa L: In Phytocannabinoids: Unraveling the Complex Chemistry and Pharmacology of Cannabis sativa. A. Douglas Kinghorn, Heinz Falk, Simon Gibbons, Jun'ichi Kobayashi. (Eds). Springer International Publishing: Cham, Switzerland. 2017, 103, 1-36.
29. M Cho, I So, JN Chun, JH Jeon, The Antitumor Effects of Geraniol: Modulation of Cancer Hallmark Pathways (Review). International Journal of Oncology. 2016, 48(5), 1772–1782. <https://doi.org/10.3892/ijo.2016.3427> <<https://doi.org/10.3892/ijo.2016.3427>> .
30. D Stevens, CJ Hudalla, M Twohig, K Organtini. Terpenes in Hemp and Cannabis Determined Using

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<https://www.waters.com/nextgen/us/en/library/application-notes/2021/analysis-of-cannabinoids-in-cannabis-plant-materials-and-edible-products-using-ultraperformance-liquid-chromatography-uplc-with-pda-and-mass-detection.html>> .

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