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1. Brickman AM, Khan UA, Provenzano FA, Yeung LK, Suzuki W, Schroeter H, Wall M, Sloan RP, Small SA. Enhancing Dentate Gyrus Function with Dietary Flavanols Improves Cognition in Older Adults. *Nat. Neurosci.*, 2014, 17(12): 1798–803.
2. Heiss C, Sansone R, Karimi H, Krabbe M, Schuler D, Rodriguez-Mateos A, Kraemer T, Cortese-Krott MM, Kuhnle GG, Spencer JP, Schroeter H, Merx MW, Kelm M, Flaviola Consortium, E.U.t.F.P. Impact of Cocoa Flavanol Intake On Age-Dependent Vascular Stiffness In Healthy Men: A Randomized, Controlled, Double-Masked Trial. *Age (Dordr)*, 2015, 37(3): 9794.
3. Kelm MA, Johnson JC, Robbins RJ, Hammerstone JF, Schmitz HH. High-Performance Liquid Chromatography Separation and Purification of Cacao (*Theobroma cacao* L.) Procyanidins According to Degree of Polymerization Using a Diol Stationary Phase. *J. Agric. Food Chem.*, 2006, 54(5): 1571–6.
4. Adamson GE, Lazarus SA, Mitchell AE, Prior RL, Cao G, Jacobs PH, Kremers BG, Hammerstone JF, Rucker RB, Ritter KA, Schmitz HH. HPLC Method for The Quantification of Procyanidins in Cocoa and Chocolate Samples and Correlation to Total Antioxidant Capacity. *J. Agric. Food Chem.*, 1999, 47(10): 4184–4188.
5. Hammerstone JF, Lazarus SA, Mitchell AE, Rucker R, Schmitz HH. Identification of Procyanidins in Cocoa (*Theobroma cacao*) and Chocolate Using High-Performance Liquid Chromatography/Mass Spectrometry. *J. Agric. Food Chem.*, 1999, 47(2): 490–496.
6. Lazarus SA, Adamson GE, Hammerstone JF, Schmitz HH. High-Performance Liquid Chromatography/Mass Spectrometry Analysis of Proanthocyanidins in Foods and Beverages. *J. Agric. Food Chem.*, 1999, 47(9): 3693–701.
7. AOAC. 2020: Official Method of Analysis 2020.05.
8. Bussy U, May BR, Olanrewaju Y, Hewitt G, Anderson N, Crozier A, Ottaviani JI, Kwik-Urbe C. Reliable, Accessible, and Transferable Method for the Quantification of Flavanols and Procyanidins in Foodstuffs and Dietary Supplements. *Food Funct.*, 2020, 11(1): 131–138.
9. Bussy U, Hewitt G, Olanrewaju Y, May BR, Anderson N, Ottaviani J, Kwik-Urbe C. Single-Laboratory Validation for the Determination of Cocoa Flavanols and Procyanidins (by Degree of

Polymerization DP1–7) in Cocoa-Based Products by Hydrophilic Interaction Chromatography Coupled with Fluorescence Detection: First Action 2020.05. *J. AOAC Int.*, 2020.

10. McCalley DV. A Study of Column Equilibration Time in Hydrophilic Interaction Chromatography. *J. Chromatogr. A*, 2018, 1554: 61–70.
11. Rimmer CA, Lippa KA, Yen J, Bussy U, Anderson N, Kwik-Urbe C. Reference Material 8403 Cocoa Flavanol Extract. *NIST Special Publication*, 2020.
12. Bell DS, Shollenberger DL. Investigation of Re-equilibration in Hydrophilic Interaction Liquid Chromatography. *LCGC*, 2016 29(12): 687–692.

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Ugo Bussy (Mars, Incorporated)、Jinchuan Yang (沃特世公司)、Hong You (Eurofins Scientific)、Nicholas Anderson (Mars, Incorporated)、Emily R Britton (沃特世公司) 以及 Catherine Kwik-Urbe (Mars, Incorporated)。

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