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1. Parr MK and Joseph JF. NDMA Impurity in Valsartan and other Pharmaceutical Products: Analytical Methods for the Determination of N-nitrosamines. JPBA 2018;164(2019):536–549.
2. US Food and Drug Administration. Liquid Chromatography-High Resolution Mass Spectrometry (LC-HRMS) Method for the Determination of Six Nitrosamine Impurities in ARB Drugs. 21 May 2019.
3. Sörgel F et al. The Contamination of Valsartan and other Sartans, Part 1: New findings. JPBA 2019;172(2019):395-405.
4. <https://www.ema.europa.eu/en/human-regulatory/post-authorisation/referral-procedures/nitrosamine-impurities> <<https://www.ema.europa.eu/en/human-regulatory/post-authorisation/referral-procedures/nitrosamine-impurities>> .Accessed 13 July 2020.
5. White CM. Understanding and Preventing (N-Nitrosodimethylamine) NDMA Contaminating of Medications. AOP 2019;54(6):611–614.
6. <https://www.fda.gov/drugs/drug-safety-and-availability/fda-updates-and-press-announcements-ndma-metformin> <<https://www.fda.gov/drugs/drug-safety-and-availability/fda-updates-and-press-announcements-ndma-metformin>> .Accessed 13 July 2020.
7. US Food and Drug Administration. Liquid Chromatography-High Resolution Mass Spectrometry (LC-HRMS) Method for the Determination of NDMA in Ranitidine Drug Substance and Drug Product. 13 September 2019.
8. Waters Corporation. XS Collision Cell White Paper. Waters Corporation White Paper [720005071EN](https://www.waters.com/waters/library.htm?cid=511436&lid=134800563&lcid=134800562&locale=en_US) <https://www.waters.com/waters/library.htm?cid=511436&lid=134800563&lcid=134800562&locale=en_US> .2014 June.
9. Tomczyk N et al. Targeted High Resolution Quantification with ToF-MRM and HD-MRM. Waters Application Note [720004728EN](https://www.waters.com/webassets/cms/library/docs/720004728en.pdf) <<https://www.waters.com/webassets/cms/library/docs/720004728en.pdf>> .2013 June.
10. US Food and Drug Administration. Control of Nitrosamine Impurities in Human Drugs: Guidance for Industry. September 2020.
11. Yang J et al. A Cautionary Tale: Quantitative LC-HRMS Analytical Procedures for the Analysis of N-Nitrosodimethylamine in Metformin. The AAPS Journal 22:89 (2020).

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