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---

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1. Horvath CG, Preiss BA, Lipsky SR. Fast Liquid Chromatography: An Investigation of Operating Parameters and the Separation of Nucleotides on Pellicular Ion Exchangers. *Analytical Chemistry*. (39) 1967, 1422.
2. Manns D, Mansfield AK. A Core-Shell Column Approach to a Comprehensive High-Performance Liquid Chromatography Phenolic Analysis of *Vitis vinifera* L. and Interspecies Hybrid Grape Juices, Wines, and Other Matrices Following either Solid Phase Extraction or Direct Injection. *Journal of Chromatography A* .(1251) 2012.111–121.
3. Nguyen J, Swann T. USP Method Modernization Using “Equivalent  $L/dp$ ” and “Equivalent  $N$ ” Allowed Changes with CORTECS  $C_8$  and CORTECS UPLC  $C_8$  Columns. Waters Application Note. 720005666EN <<https://www.waters.com/nextgen/us/en/library/application-notes/2016/usp-method-modernization-using-equivalent-l-dp-and-equivalent-n-allowed-changes.html>> , 2016.
4. Berthelette K, Swann T. 使用CORTECS  $C_8$  色谱柱根据USP方法分析氯雷他定中的杂质. 沃特世应用纪要. 720005579EN <<https://www.waters.com/nextgen/us/en/library/application-notes/2016/usp-method-for-the-impurity-analysis-of-loratidine-using-cortecs-c8-columns.html>> , 2016.
5. Nguyen J, Swann T. Examples of USP Method Modernizations Using “Equivalent  $L/dp$ ” and “Equivalent  $N$ ” Allowed Changes with CORTECS  $C_8$  Columns. Waters Application Note. 720005988EN <<https://www.waters.com/nextgen/us/en/library/application-notes/2017/usp-using-equivalent-l-dp-and-equivalent-n-allowed-changes-with-cortecs-c8-columns.html>> , 2017.
6. Cabooter D, Fanigliulo A, Bellazzi G, Allieri B, Rottigni A, Desmet G. Relationship Between the Particle Size Distribution of Commercial Fully Porous and Superficially Porous High-Performance Liquid Chromatography Column Packings and Their Chromatographic Performance. *Journal of Chromatography A*.(1217) 2010.7074-7081.
7. DeStefano J, Langlois T, Kirkland J. Characteristics of Superficially Porous Silica Particles for Fast HPLC: Some Performance Comparisons with Sub-2- $\mu\text{m}$  Particles. *Journal of Chromatographic Science* (46) 2008.254-260.
8. Waters T, Wyndham K, Shiner S, Walsh D, Okandeji O, Savaria M, Iraneta P, Turner J, Boissel C, McCall

S. [Advances in Particle and Packing Technology for Higher Efficiency Ultraperformance LC Columns <https://www.waters.com/webassets/cms/library/docs/2013hplc\\_walter\\_packing\\_technology.pdf>](https://www.waters.com/webassets/cms/library/docs/2013hplc_walter_packing_technology.pdf)  
.Waters Scientific Poster HPLC Conference 2013.

9. Hayes R, Ahmed A, Edge T, Zhang H. Core-Shell Particles: Preparation, Fundamentals, and Applications in High Performance Liquid Chromatography. *Journal of Chromatography A*.(1357) 2014.36-52.

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Jennifer Simeone, Pat McConville, Mark Trahan

---

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[Empower色谱数据系统 <https://www.waters.com/10190669>](https://www.waters.com/10190669)

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© 2024 Waters Corporation. All Rights Reserved.

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