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1. Wolever Th, A. Piekarz, M. Hollands, K. Younker, Rapaille.A, Goosens.J, etc.2002. Sugar Alcohols and Diabetes: A Review.*Canadian Journal of Diabetes*.26(4): 356–362.
2. Rapaille.A, Goosens.J, & Heume.M. (2016) Sugar alcohols.*Encyclopedia of Food and Health*, 211–216.
3. Ruskone-Fourmestraux.A, Attar.A, Chassard.D, Coffin.B, Bornet.F, Bouhnik.Y. (2003) A Digestive Tolerance Study of Maltitol after Occasional and Regular Consumption in Healthy Humans.*European Journal of Clinical Nutrition*, 57, 26–30.
4. Akram Hossain; Fuminori Yamaguchi; Tatsuhiro Matsuo; Ikuko Tsukamoto; Yukiyasu Toyoda; Masahiro Ogawa; Yasuo Nagata; Masaaki Tokuda.(November 2015) "Rare Sugar d-allulose: Potential Role and Therapeutic Monitoring in Maintaining Obesity and Type 2 Diabetes Mellitus".*Pharmacology & Therapeutics*.155: 49–59.
5. U.S. FDA (2020) Guidance for Industry: The Declaration of Allulose and Calories from Allulose on Nutrition and Supplement Facts Labels.<https://www.fda.gov/regulatory-information/search-fda-guidance-documents/guidance-industry-declaration-allulose-and-calories-allulose-nutrition-and-supplement-facts-labels> <<https://www.fda.gov/regulatory-information/search-fda-guidance-documents/guidance-industry-declaration-allulose-and-calories-allulose-nutrition-and-supplement-facts-labels>> .
6. Southey, F. (2021) Allulose approval in Europe to be sought by New Ingredients Consortium.*Food Navigator*.<https://www.foodnavigator.com/Article/2021/12/07/Allulose-approval-in-Europe-to-be-sought-by-new-ingredients-consortium> <<https://www.foodnavigator.com/Article/2021/12/07/Allulose-approval-in-Europe-to-be-sought-by-new-ingredients-consortium>> .

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Arc HPLC系统 <<https://www.waters.com/waters/nav.htm?cid=135068659>>

[2414示差折光\(RI\)检测器 <https://www.waters.com/waters/514425>](https://www.waters.com/waters/514425)

[Empower色谱数据系统 <https://www.waters.com/10190669>](https://www.waters.com/10190669)

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© 2022 Waters Corporation. All Rights Reserved.

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