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Christian Lanshoeft, Frederic Lozac'h - Novartis Pharma AG。

---

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1. ICH M3 (R2) [Non-clinical Safety Studies for the Conduct of Human Clinical Trials for Pharmaceuticals](https://www.ema.europa.eu/en/ich-m3-r2-non-clinical-safety-studies-conduct-human-clinical-trials-pharmaceuticals#current-effective-version-section) <<https://www.ema.europa.eu/en/ich-m3-r2-non-clinical-safety-studies-conduct-human-clinical-trials-pharmaceuticals#current-effective-version-section>> , 2013.
2. FDA Guidance for Industry. [Safety Testing of Drug Metabolites](https://www.fda.gov/regulatory-information/search-fda-guidance-documents/safety-testing-drug-metabolites) <<https://www.fda.gov/regulatory-information/search-fda-guidance-documents/safety-testing-drug-metabolites>> , 2020.
3. Broeckling C, Yao L, Isaac G, Gioioso M, Ianchis V, Vissers JPC. Application of Predicted Collisional Cross Section to Metabolome Databases to Probabilistically Describe the Current and Future Ion Mobility Mass Spectrometry. *J Am Soc Mass Spectrom*, 2021, 32, 661–669.
4. Higton D, Palmer ME, Vissers JPC, Mullin LG, Plumb RS, Wilson ID. Use of Cyclic Ion Mobility Spectrometry (cIM)-Mass Spectrometry to Study the Intramolecular Transacylation of Diclofenac Acyl Glucuronide. *Anal Chem*, 2021, 93, 20, 7413–7421.
5. Connolly JFRB, Munoz-Muriedas J, Laphorn C, Higton D, Vissers JPC, Webb A, Beaumont C, Dear G.J. Investigation into Small Molecule Isomeric Glucuronide Metabolite Differentiation Using In Silico and Experimental Collision Cross-Section Values. *J. Am.Soc.Mass Spectrom.* 2021; 32 (8): 1976–1986.
6. Holdsworth C, Clayton R, Robinson H, Lord-Mears C, Kendrick J. Utilization of Ion Mobility Enabled

Collisional Cross Section Measurements for the Comparison of Metabolites across Differing Chromatographic Methods. Poster DMDG 2016.

7. Giles K, Ujma J, Wildgoose J, Pringle S, Richardson K, Langridge D, Green M. A Cyclic Ion Mobility-Mass Spectrometry System. *Anal. Chem.* 2019, 91, 13, 8564–8573.

---

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[ACQUITY UPLC I-Class PLUS系统 <https://www.waters.com/134613317>](https://www.waters.com/134613317)

[SELECT SERIES Cyclic IMS <https://www.waters.com/waters/nav.htm?cid=135021297>](https://www.waters.com/waters/nav.htm?cid=135021297)

[SYNAPT G2-Si质谱仪 <https://www.waters.com/134740653>](https://www.waters.com/134740653)

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[MassLynx MS软件 <https://www.waters.com/513662>](https://www.waters.com/513662)

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© 2022 Waters Corporation. All Rights Reserved.

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