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1. Uddin SMK *et al.* Halal and Kosher Gelatin: Applications as Well as Detection Approaches with Challenges and Prospects. *Food Bioscience* (2021) 41:101422.
2. Ali ME, Suleiman N. Eleven Shades of Food Integrity: A Halal Supply Chain Perspective. *Trends in Food*

Sci. Technol.(2018) 71:216–224.

3. Hameed AM *et al.*A Review of Gelatin Source Authentication Methods.*Trop.Life Sci.Res.*(2018), 29(2):213–227.
4. Rohman A *et al.*Review on Analytical Methods for Analysis of Porcine Gelatine in Food and Pharmaceutical Products for Halal Authentication.*Trends in Food Sci.Technol.* (2020) 101:122–132.
5. Ishaq A *et al.*Potentiality of Analytical Approaches to Determine Gelatin Authenticity in Food Systems: A Review.*LWT* (2020) 121:108968.
6. Valetta M *et al.*Mass Spectrometry-Based Protein and Peptide Profiling for Food Frauds, Traceability, and Authenticity Assessment. *Food Chem.* (2021) 365:130456.
7. Cheng X-L *et al.*Identification of Five Gelatins by Ultra-Performance Liquid Chromatography/Time-Of-Flight Mass Spectrometry (UPLC/QToF-MS) Using Principal Component Analysis.*J. Pharm.Biomed.Anal.* (2012) 62:191-195.
8. Grundy H *et al.* A Mass Spectrometry Method for the Determination of the Species of Origin of Gelatine in Foods and Pharmaceutical Products.*Food Chem.* (2016) 190:276-284.
9. Sha M-X *et al.*The Identification of Three Mammalian Gelatins by Liquid Chromatography-High Resolution Mass Spectrometries.*LWT* (2018) 89:74e86.
10. Deng G *et al.*Recent Advances in Animal Origin Identification of Gelatin-based Products Using Liquid Chromatography-Mass Spectrometry Methods: A Mini Review. *Reviews in Anal.Chem.*(2020) 9:260-271.

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© 2023 Waters Corporation. All Rights Reserved.