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---

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1. M. Anastassiades, S. J. Lehotay, D. Stajnbaher, F. J. Schenck, Fast and Easy Multiresidue Method Employing Acetonitrile Extraction/Partitioning and “Dispersive Solid-Phase Extraction” for the Determination of Pesticide Residues in Produce, *J. AOAC Int.*, 86 (2003) 412–431.
2. Dutch mini-Luke ( “NL-” ) Extraction Method Followed by LC and GC-MS/MS for Multiresidue Analysis of Pesticides in Fruits and Vegetables, URL: <https://www.eurl-pesticides.eu/userfiles/file/NL-miniLuke-extraction-method.pdf> <<https://www.eurl-pesticides.eu/userfiles/file/NL-miniLuke-extraction-method.pdf>> .
3. M. Anastassiades; A.-K. Wachtler; D. I. Kolberg; E. Eichhorn; H. Marks; A. Benkenstein; S. Zechmann; D. Mack; C. Wi Idgrube; A. Barth; I. Sigalov; S. Görlich; D. Dörk and G. Cerchia. Quick Method for the Analysis of Highly Polar Pesticides in Food Involving Extraction with Acidified Methanol and LC - or IC MS/MS Measurement - I. *Food of Plant Origin* (QuPPE-PO-Method) – 第12版 (2021年7月23日在EURL-SRM网站上发布) ; URL: [https://www.eurl-pesticides.eu/docs/public/tmpl\\_article.asp?CntID=887&LabID=200&Lang=EN](https://www.eurl-pesticides.eu/docs/public/tmpl_article.asp?CntID=887&LabID=200&Lang=EN) <[https://www.eurl-pesticides.eu/docs/public/tmpl\\_article.asp?CntID=887&LabID=200&Lang=EN](https://www.eurl-pesticides.eu/docs/public/tmpl_article.asp?CntID=887&LabID=200&Lang=EN)> .
4. Hird S, Adams S, De-Alwis J. Evaluation of the Performance of a Method for the Determination of Highly Polar, Anionic Pesticides in Foodstuffs Using LC-MS/MS, Waters Application Note, [720007505](#), 2022.
5. De-Alwis J, Williams J, Hird S, Adams S. 使用实验室间研究评估LC-MS/MS方法测定农作物和食品中阴离子型极性农药残留的性能, 沃特世应用纪要, [720007154ZH](#), 2021.
6. Kumar K P, Bhaska K, Gorella T, Wagh P. 使用UPLC-MS和阴离子型极性农药分析色谱柱测定葡萄中的阴离子型极性农药. 沃特世应用纪要, [720006925ZH](#), 2020.
7. Ross E, De-Alwis J, Adams S, Williams J, Shah D D. Determination of Anionic Polar Pesticides in High Water Foodstuffs. Waters Application Note, [720006645](#), 2019.

---

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[Xevo TQ Absolute </nextgen/cn/zh/products/mass-spectrometry-systems/xevo-tq-absolute.html>](/nextgen/cn/zh/products/mass-spectrometry-systems/xevo-tq-absolute.html)

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© 2022 Waters Corporation. All Rights Reserved.

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