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1. Niu Y *et al.* Atmospheric Pressure Chemical Ionization Source as an Advantageous Technique for Gas Chromatography-Tandem Mass Spectrometry. *Trends Anal.Chem.* (2020) 132:116053.
2. Cherta L *et al.* Application of Gas Chromatography-(Triple Quadrupole) Mass Spectrometry With Atmospheric Pressure Chemical Ionization for the Determination of Multiclass Pesticides in Fruits and

Vegetables. *J Chromatogr.A* (2013) 1314:224–240.

3. Saito-Shida S *et al.* Quantitative Analysis of Pesticide Residues in Tea by Gas Chromatography–Tandem Mass Spectrometry With Atmospheric Pressure Chemical Ionization. *J Chromatogr.B* (2020a) 1143:122057.
4. Saito-Shida S *et al.* Multi-Residue Determination of Pesticides in Green Tea by Gas Chromatography-Tandem Mass Spectrometry With Atmospheric Pressure Chemical Ionisation Using Nitrogen as the Carrier Gas. *Food Addit.Contam.Part A* (2020b) 38(1): 125–135.
5. European Committee for Standardisation (CEN) EN 15662:2018. Foods of Plant Origin - Multimethod for the Determination of Pesticide Residues Using Gc- And LC- Based Analysis Following Acetonitrile Extraction/Partitioning and Clean-up by Dispersive Spe - Modular Quechers-Method.
6. Document No.SANTE/12682/2019.Guidance Document on Analytical Quality, Control, and Method Validation Procedures for Pesticides Residues Analysis in Food and Feed.2019.
7. Fussell R *et al.* Assessment of the Stability of Pesticides during Cryogenic Sample Processing.1. Apples. *Agric.Food Chem.*(2002) 50(3):441–448.
8. Rutkowska E *et al.* Modification of Multiresidue QuEChERS Protocol to Minimize Matrix Effect and Improve Recoveries for Determination of Pesticide Residues in Dried Herbs Followed by GC-MS/MS. *Food Anal.Methods* (2018) 11:709–724.

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MassLynx MS软件 <<https://www.waters.com/513662>>

TargetLynx <<https://www.waters.com/513791>>

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© 2022 Waters Corporation. All Rights Reserved.

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