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1. Guo S, Xu X, Zhou X and Huang Y. A rapid and simple UPLC-MS/MS method using collagen marker peptides for identification of porcine gelatin. *RSC Adv.* (2018), 8: 3768–3773.
2. Shabani H, Mehdizadeh M, Mousavi S M, Dezfouli E A, Solgi T, Khodaverdi M, Rabiei M, Rastegar H and Alebouyeh M. Halal authenticity of gelatin using species-specific PCR. *Food Chem.* (2015), 184: 203-206.
3. Uddin SMK et al. Halal and Kosher Gelatin: Applications as Well as Detection Approaches with Challenges and Prospects. *Food Bioscience* (2021) 41:101422.
4. Nhari R M H R, Ismail A, and Man Y B C. Analytical Methods for Gelatin Differentiation from Bovine and Porcine Origins and Food Products. *J. Food Sci.* (2012), 71(1): R42–46.
5. Deng G *et al.* Recent Advances in Animal Origin Identification of Gelatin-Based Products Using Liquid Chromatography-Mass Spectrometry Methods: A mini review. *Reviews in Analytical Chemistry* (2020) 39:260–271.
6. Chantel Lee, Li Yan Chan, A Complete Discovery Workflow for Species-Specific Gelatin Identification, Waters Application Note, [720007533](https://www.waters.com/720007533), 2022.
7. Song E, Gao Y, Wu C, Shi T, Nie S, Fillmore T L, Schepmoes A A, Gritsenko M A, Qian W, Smith R D, Rodland K D, Liu T. Targeted Proteomic Assays For Quantitation of Proteins Identified by Proteogenomic Analysis of Ovarian Cancer. *Sci.Data* (2017), 1–13.

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[MassLynx MS软件 <https://www.waters.com/513662>](https://www.waters.com/513662)

[TargetLynx <https://www.waters.com/513791>](https://www.waters.com/513791)

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© 2022 Waters Corporation. All Rights Reserved.

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