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1. Schrenk D, Bignami M, *et al.* EFSA Panel on Contaminants in the Food Chain (EFSA CONTAM Panel), Risk to Human Health Related to the Presence of Perfluoroalkyl Substances in Food. *EFSA J.* 2020;18(9).
2. Treat SA. Forever Chemicals and Agriculture Case Study: Maine Accelerates Across-The-Board Action to Address Pfas Chemicals Harming Farmers and Rural Communities. Institute for Agriculture and Trade Policy. November 2021.
3. Maine Department of Agriculture, Conservation and Forestry.
<https://www.maine.gov/dacf/ag/pfas/index.shtml#impact> <
<https://www.maine.gov/dacf/ag/pfas/index.shtml#impact>> . Accessed 15 June 2022.
4. Sznajder-Katarzynska K, Surma M, Wiczowski W, Cieslik E. The Perfluoroalkyl Substance (Pfas) Contamination Level in Milk and Milk Products in Poland. *International Dairy Journal.* 96, 2019, 73-84.
5. Macheke LR, Olowoyo JO, Mugivhisa LL, Afafe OA. Determination and assessment of human dietary intake of per and polyfluoroalkyl substances in retail dairy milk and infant formula from South Africa. *Science of the Total Environment.* 755, 2021.
6. Barbarossa A, Gazzotti T, Zironi E, Serraino A, Pagliuca G. Monitoring the presence of perfluoroalkyl substances in Italian cow milk. *J. Dairy Sci.* 97, 2014, 3339-3343.
7. Organtini K, Hird S, Adams S. 利用QuEChERS从可食用农产品中提取全氟烷基和多氟烷基化合物(PFAS)并在Xevo TQ-XS上进行高灵敏度分析. 沃特世应用纪要 [720007333ZH](#). 2021年8月.

8. Food and Drug Administration Foods Program.Guidelines for the Validation of Chemical Methods in Food, Feed, Cosmetics, and Veterinary Products, 3rd edition.October 2019.
9. European Union Reference Laboratory.Guidance Document on Analytical Parameters for the Determination of Per- and Polyfluoroalkyl Substances (PFAS) in Food and Feed.Version 1.2. May 2022.

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© 2022 Waters Corporation. All Rights Reserved.

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[Cookie](#)

[Cookie](#)

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