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1. Moser J., Schlingemann J., Saal C. N-Nitrosamines Impurities in Pharmaceuticals the Abrupt Challenges That Resulted, the Evolving Science, and the Regulatory Framework. *Journal of Pharmaceutical Sciences* 2023, 112 (5), 1161–1162. DOI: <https://doi.org/10.1016/j.xphs.2023.01.016> <<https://doi.org/10.1016/j.xphs.2023.01.016>> .
2. Teasdale A., Moser J., Ford J. G., Creasy, J. N -Nitrosamines. In *Mutagenic Impurities*, Teasdale, A. Ed., Wiley, 2021; pp 269–320.
3. Bharate, S. S. Critical Analysis of Drug Product Recalls Due to Nitrosamine Impurities. *Journal of Medicinal Chemistry* 2021, 64 (6), 2923–2936. DOI: 10.1021/acs.jmedchem.0c02120.
4. Lame M. E., Hatch, L. High Sensitivity Quantitation of Nitrosamine Genotoxic Impurities: LC-MS Analysis of Ranitidine Drug Product using the Waters ACQUITY UPLC I-Class/Xevo TQ-XS Tandem Quadrupole Mass Spectrometer. Waters Application Note. [720006899](https://www.waters.com/content/dam/waters/technical-content/application-notes/720006899.pdf), 2020.
5. Maziarz M., Rainville P. D., Bartlett A. 利用Xevo™ TQ Absolute串联四极杆质谱仪和ACQUITY™ Premier系统对二甲双胍中的亚硝酸胺进行高灵敏度定量分析, 沃特世应用纪要, [720007725ZH](https://www.waters.com/content/dam/waters/technical-content/application-notes/720007725ZH.pdf), 2022.
6. Chang S.-H., Ho H.-Y., Chang C.-C., Zang C.-Z., Hsu Y.-H., Lin M.-C., Tseng S.-H., Wang D.-Y. Evaluation and Optimization of a HS-SPME-assisted GC-MS/MS Method for Monitoring Nitrosamine Impurities in Diverse Pharmaceuticals. *Journal of Pharmaceutical and Biomedical Analysis* 2022, 221, 115003. DOI: <https://doi.org/10.1016/j.jpba.2022.115003> <<https://doi.org/10.1016/j.jpba.2022.115003>> .
7. Wichitnithad W., Nantaphol S., Noppakhunsomboon K., Thitikornpong W., Rojsitthisak, P. Current Status and Prospects of Development of Analytical Methods for Determining Nitrosamine and N-Nitroso Impurities in Pharmaceuticals. *Talanta* 2023, 254, 124102. DOI: <https://doi.org/10.1016/j.talanta.2022.124102> <<https://doi.org/10.1016/j.talanta.2022.124102>> .
8. David F., Vanhoenacker G., Sandra K., Sandra P., Bristow T., Harrison M. Strategic Approaches to the Chromatographic Analysis of Mutagenic Impurities. In *Mutagenic Impurities*, Teasdale, A. Ed., Wiley, 2021; pp 381–438.
9. Zhao Y.-Y., Boyd J., Hrudey S. E., Li X.-F. Characterization of New Nitrosamines in Drinking Water Using Liquid Chromatography Tandem Mass Spectrometry. *Environmental Science & Technology* 2006, 40 (24), 7636–7641. DOI: 10.1021/es061332s.
10. Nagendla N. K., Shaik H., Bala Subrahanyam S., Godugu D., Reddy Mudiam, M. K. Development,

Validation, and Estimation of Measurement Uncertainty for the Quantitative Determination of Nitrosamines in Sartan Drugs Using Liquid Chromatography-Atmospheric Pressure Chemical Ionization-Tandem Mass Spectrometry. *Journal of Chromatography Open* 2022, 2, 100053. DOI: <https://doi.org/10.1016/j.jcoa.2022.100053> <<https://doi.org/10.1016/j.jcoa.2022.100053>> .

11. Lee J.-H., Lee S.-U., Oh J.-E. Analysis of Nine Nitrosamines in Water by Combining Automated Solid-Phase Extraction With High-Performance Liquid Chromatography-Atmospheric Pressure Chemical Ionisation Tandem Mass Spectrometry. *International Journal of Environmental Analytical Chemistry* 2013, 93 (12), 1261–1273. DOI: 10.1080/03067319.2013.803281.
12. Ripollés C., Pitarch E., Sancho J. V., López F. J., Hernández F. Determination of Eight Nitrosamines in Water at the ng L⁻¹ Levels by Liquid Chromatography Coupled to Atmospheric Pressure Chemical Ionization Tandem Mass Spectrometry. *Analytica Chimica Acta* 2011, 702 (1), 62–71. DOI: <https://doi.org/10.1016/j.aca.2011.06.024>.

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ACQUITY Premier系统 <<https://www.waters.com/waters/nav.htm?cid=135077739>>

Xevo TQ Absolute <<https://www.waters.com/nextgen/global/products/mass-spectrometry-systems/xevo-tq-absolute.html>>

MassLynx MS软件 <<https://www.waters.com/513662>>

TargetLynx <<https://www.waters.com/513791>>

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© 2023 Waters Corporation. All Rights Reserved.

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