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1. Ly T-K *et al.* Quantification of 397 Pesticide Residues in Different Types of Commercial Teas: Validation of High Accuracy Methods and Quality Assessment. *Food Chem.* (2022) 370:130986.
2. Comprehensive Strategy for Pesticide Residue Analysis in Cocoa Beans Through Qualitative and Quantitative Approach. *Food Chem.* (2022) 368:130778.
3. European Committee for Standardisation (CEN) EN 15662:2018. Foods of Plant Origin - Multimethod for the Determination of Pesticide Residues Using GC- And LC- Based Analysis Following Acetonitrile Extraction/Partitioning and Clean-up by Dispersive SPE - Modular QuEChERS-Method.
4. Shah D., Wood J., Fujimoto G., McCall E., Hird S., Hancock P. 使用UPLC-MS/MS定量测定水果、蔬菜、谷类和红茶中农药的多残留分析方法. 沃特世应用纪要. [720006886ZH](#). 2020.
5. Young M., Shia J. Oasis PRiME HLB Cartridges Now Available in Syringe Compatible Plus Format. Waters Application Note. [720006017](#). 2017.
6. Document No. SANTE/11312/2021. [Guidance Document on Analytical Quality, Control, and Method Validation Procedures for Pesticides Residues Analysis in Food and Feed. 2021](#) <https://www.eurl-pesticides.eu/userfiles/file/EurlALL/SANTE_11312_2021.pdf> .

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[ACQUITY UPLC I-Class PLUS系统 <https://www.waters.com/134613317>](https://www.waters.com/134613317)

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MassLynx MS软件 <<https://www.waters.com/513662>>

TargetLynx <<https://www.waters.com/513791>>

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© 2024 Waters Corporation. All Rights Reserved.

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