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1. Significance of Lipolysis in the Manufacture and Storage of Dairy Products. IDF Bulletin No.144. Brussels, Belgium: *International Dairy Federation*. IDF (1987).
2. Clarke, H. J., McCarthy, W. P., O'Sullivan, M. G., Kerry, J. P., and Kilcawley, K. N., Oxidative Quality of Dairy Powders: Influencing Factors and Analysis. *Foods (Basel, Switzerland)*, 2021. 10, 2315. <https://doi.org/10.3390/foods10102315> <<https://doi.org/10.3390/foods10102315>> .
3. Imamura, F., Fretts, A., Marklund, M., Ardisson Korat, A. V., Yang, W. S., Lankinen, M., Qureshi, W., Helmer, C., Chen, T. A., Wong, K., Bassett, J. K., Murphy, R., Tintle, N., Yu, C. I., Brouwer, I. A., Chien, K. L., Frazier-Wood, A. C., Del Gobbo, L. C., Djoussé, L., Geleijnse, J. M., Fatty Acids and Outcomes Research Consortium (FORCE). Fatty Acid Biomarkers of Dairy Fat Consumption and Incidence of Type 2 Diabetes: A Pooled Analysis of Prospective Cohort Studies. *PLoS medicine*, 15, e1002670. 2018. <https://doi.org/10.1371/journal.pmed.1002670> <<https://doi.org/10.1371/journal.pmed.1002670>> .
4. de Oliveira Otto, M. C., Lemaitre, R. N., Song, X., King, I. B., Siscovick, D. S., Mozaffarian, D., 2018. Serial Measures of Circulating Biomarkers of Dairy Fat and Total and Cause-Specific Mortality in Older Adults: the Cardiovascular Health Study. *The American journal of clinical nutrition*, 108, 476–484.

<https://doi.org/10.1093/ajcn/nqy117> <<https://doi.org/10.1093/ajcn/nqy117>> .

5. Pilar Gómez-Cortés, Manuela Juárez, Miguel Angel de la Fuente, 2018. Milk Fatty Acids and Potential Health Benefits: An Updated Vision. *Trends in Food Science & Technology*, 81, 1–9.

<https://doi.org/10.1016/j.tifs.2018.08.014> <<https://doi.org/10.1016/j.tifs.2018.08.014>> .

6. Kokotou, M. G., Mantzourani, C., & Kokotos, G., Development of a Liquid Chromatography-High Resolution Mass Spectrometry Method for the Determination of Free Fatty Acids in Milk. *Molecules* (Basel, Switzerland). 2020, 25, 1548. <https://doi.org/10.3390/molecules25071548> <

<https://doi.org/10.3390/molecules25071548>> .

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ACQUITY UPLC I-Class PLUS系统 <<https://www.waters.com/134613317>>

ACQUITY QDa II质谱检测器 <<https://www.waters.com/nextgen/global/products/mass-spectrometry/mass-spectrometry-systems/acquity-qda-ii-mass-detector.html>>

MassLynx质谱软件 <<https://www.waters.com/513662>>

TargetLynx <<https://www.waters.com/513791>>

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© 2024 Waters Corporation. All Rights Reserved.

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