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1. Schymanski E, Zhang J, Thiessen P, Chirsir P, Kondic, T, Bolton, E. Per- and Polyfluoroalkyl substances (PFAS) in PubChem: 7 Million and Growing. *Environ.Sci.Technol.* 2023 57 (44): 16918–16928
2. US Environmental Protection Agency. *CompTox Chemicals Dashboard*. PFAS Structures in DSSTox (2022年8月更新) .<https://comptox.epa.gov/dashboard/chemical-lists/PFASSTRUCTV5> <
<https://comptox.epa.gov/dashboard/chemical-lists/PFASSTRUCTV5>>
3. US Environmental Protection Agency. Analysis of Per- and Polyfluoroalkyl Substances (PFAS) in Aqueous, Solid, Biosolids, and Tissue Samples by LC-MS/MS. January 2024
4. ASTM D8421–22, Standard Test Method for Determination of Per- and Polyfluoroalkyl Substances (PFAS) in Aqueous Matrices by Co-solvation followed by Liquid Chromatography Tandem Mass Spectrometry (LC/MS/MS), ASTM International, West Conshohocken, PA, 2022, www.astm.org <
<http://www.astm.org/>>
5. Liu Y, D’ Agostino LA, Qu G, Jiang G, Martin JW. High-resolution mass spectrometry (HRMS) methods for non-target discovery and characterization of poly- and perfluoroalkyl substances (PFASs) in environmental and human samples. *Trends Anal Chem* 2019, 115420
6. McCullagh M, Hodgkinson M, Dilllon L, van Bevel B, Burgess JA, Jogsten IE. 使用UPLC/带有离子淌度的质谱检测和鉴定环境样品中的全氟烷基和全氟烷基化合物(PFAS). 沃特世应用纪要, [720005090ZH](https://www.waters.com). 2014年6月.
7. Dodds JN, Hopkins ZR, Knappe DRU, Baker ES. Rapid Characterization of Per- and Polyfluoroalkyl substances (PFAS) by Ion Mobility Spectrometry-Mass Spectrometry (IMS-MS). *Anal Chem* 2020, 92, 6, 4427–4435

8. Organtini KL, Rosnack KJ, Plummer C, Hancock P, Burt O. 根据EPA 1633分析全氟烷基和多氟烷基化合物(PFAS)第2部分：水相基质分析.沃特世应用纪要, [720008143ZH](#).2023年12月.
9. Organtini KL, Rosnack KJ, Hancock P. 使用Atlantis BEH C18 AX混合模式色谱柱将超短链PFAS纳入单次进样的PFAS检测范围.沃特世应用纪要, [720008034ZH](#).2023年8月.
10. Twohig M, Fujimoto G, Mohan A, Organtini K, Rosnack KJ, Hird S. 环境样品中全氟烷基和多氟烷基化合物(PFAS)的非靶向分析方法.沃特世应用纪要 [720007184ZH](#).2021年3月.
11. Charbonnet JA, McDonough CA, Xiao F, Schwichtenberg T, Cao D, Kaserzon S, Thomas KV, Dewapriya P, Place BJ, Schymanski EL, Field JA, Helbling DE, Higgins, CP.Communicating Confidence of Per- and Polyfluoroalkyl Substance Identification via High-Resolution Mass Spectrometry. *Environ.Sci.Technol.Lett.*2022, 9, 473–481
12. Mullin L, Jobst K, DiLorenzo RA, Plumb R, Reiner EJ, Yeung LWY, Jogsten IE.Liquid chromatography-ion mobility-high resolution mass spectrometry for analysis of pollutants in indoor dust: Identification and predictive capabilities.*Analytica Chimica Acta.*2020 Aug 1125: 29–40
13. Foster M, Rainey M, Watson C, Dodds JN, Kirkwood KI, Fernandez FM, Baker ES.Uncovering PFAS and Other Xenobiotics in the Dark Metabolome Using Ion Mobility, Mass Defect Analysis, and Machine Learning.*Environ.Sci.Technol.*2022, 56, 12, 9133–9143
14. MacNeil A, Li X, Amiri R, Muir DCG, Simpson A, Simpson MJ, Dorman FL, Jobst KJ.Gas Chromatography-(Cyclic) Ion Mobility Mass Spectrometry: A Novel Platform for the Discovery of Unknown Per-/Polyfluoroalkyl Substances.*Anal.Chem.*2022, 94: 11096–11103
15. Huset CA, Barlaz MA, Barofsky DF, Field JA.Quantitative Determination of Fluorochemicals in Municipal Landfill Leachates.*Chemosphere.* 2011 Mar 82, 10, 1380–1386
16. Newton S, McMahan R, Stoeckel JA, Chislock M, Lindstrom A, Strynar M. Novel Polyfluorinated Compounds Identified Using High Resolution Mass Spectrometry Downstream of Manufacturing Facilities near Decatur, Alabama.*Environ.Sci.Technol.*2017, 51, 1544–1552

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[ACQUITY UPLC I-Class PLUS系统 <https://www.waters.com/134613317>](https://www.waters.com/134613317)

[SELECT SERIES Cyclic IMS <https://www.waters.com/waters/global/SELECT-SERIES-Cyclic-IMS-ion-mobility-mass-spectrometer/nav.htm?cid=135021297>](https://www.waters.com/waters/global/SELECT-SERIES-Cyclic-IMS-ion-mobility-mass-spectrometer/nav.htm?cid=135021297)

[waters_connect <https://www.waters.com/waters/nav.htm?cid=135040165>](https://www.waters.com/waters/nav.htm?cid=135040165)

[UNIFI科学信息系统 <https://www.waters.com/waters/global/UNIFI-Scientific-Information-System/nav.htm?cid=134801648>](https://www.waters.com/waters/global/UNIFI-Scientific-Information-System/nav.htm?cid=134801648)

[MassLynx MS软件 <https://www.waters.com/513662>](https://www.waters.com/513662)

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© 2025 Waters Corporation. All Rights Reserved.

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